

# ENC DC Series Autosampler Enclosures

ENC series DC Autosampler Enclosures protect your samples, your employees, and your peace of mind.



## DC Autosampler Enclosure Features

- Save bench space with an enclosure sized to fit your autosampler (see below)
- Flat pack design for inexpensive shipping
- Shatter-proof polycarbonate construction
- Quickly evacuates harmful fumes and protects samples from contamination
- Prevents contamination and corrosion
- Easy access to samples and autosampler components through front and rear
- Chemical resistant materials
- Multiple configuration possibilities for vent and fan
- Multiple ports for tubing and power cords
- Can be used without the fan and filter to function as a convenient dust cover
- Included 4" exhaust port flange
- Optional fan maintains positive pressure
- Optional HEPA filter
- Optional rigid base for irregular lab surfaces
- Multiple models, compatible with CETAC ASX-130, ASX-260, ASX-280, ASX-520, ASX-520HS, ASX-560, EXR-8, XLR-8, and XLR-860 autosamplers and more!

## Technical Specifications

Dimensions	Height	Width	Depth	Weight
<b>ENC-560DC</b>	65 cm (26 in)	60 cm (24 in)	58 cm (23 in)	11 kg (25 lb)
<b>ENC-280DC</b>	65 cm (26 in)	43 cm (17 in)	58 cm (23 in)	8 kg (17 lb)
<b>ENC-860DC</b>	98 cm (39 in)	111 cm (44 in)	67 cm (27 in)	22 kg (49 lb)

Shipping Details	Box 1	Box 2	Base
<b>ENC-560DC</b>	28" x 28" x 1.5" at 15 lb	8" x 8" x 24" at 10 lb	28" x 25.5" x 1" at 18 lb
<b>ENC-280DC</b>	28" x 28" x 1.5" at 15 lb	8" x 8" x 17" at 7 lb	28" x 18" x 1" at 12 lb
<b>ENC-860DC</b>	33" x 33" x 1.75" at 30 lb	12" x 12" x 48" at 23 lb	48" x 30" x 1" at 35 lb

### Materials

polycarbonate frame  
silicone and polyethylene plugs  
nylon fittings

### Exhaust

4" fume hood exhaust attachment

### Optional Accessories

HEPA filter  
Positive pressure fan  
Gas flow meter (up to 5 L per minute)  
Base for irregular surfaces

### Warranty



1 year limited (by manufacturer)



## ENC-560 DC

for ASX-560 and 4-rack autosamplers

Part Number	Description
<b>ENC-99-0010</b>	ENC-560 DC Autosampler Enclosure



## ENC-280 DC


for ASX-280 and 2-rack autosamplers

Part Number	Description
<b>ENC-99-0011</b>	ENC-280 DC Autosampler Enclosure



## ENC-860 DC



for XLR-8<sub>60</sub> and 8-rack autosamplers

Part Number	Description
<b>ENC-99-0012</b>	ENC-860 DC Autosampler Enclosure



## ENC DC Series Enclosure Spares

Part Number	Description	Image
	Enclosure Fan - positive air flow for use with HEPA filter when vent is not available. Provides a flow rate of 121 CFM.	
<b>32-0099-045</b>	North America	
<b>32-0100-045</b>	Europe	
<b>32-0101-045</b>	Australia	
	HEPA Filter (Requires adapter plate)	
<b>32-0162-045</b>	Filter with adapter plate	
<b>32-0102-045</b>	Filter replacement only	

Part Number	Description	Image
	Exhaust flange	
<b>32-0288-045</b>	2.5" Diameter	
<b>32-0289-045</b>	3.0" Diameter	
<b>32-0290-045</b>	3.5" Diameter	
<b>32-0103-045</b>	4.0" Diameter	
	Autosampler Base	
<b>32-0105-045</b>	for ASX-280	
<b>32-0104-045</b>	for ASX-560	
<b>32-0106-045</b>	for XLR-8 <sub>60</sub>	